


Product Datasheet

SARS-CoV-2 (COVID-19) Envelope Antibody (biotin) (orb1239041)

Description

SARS-CoV-2 (COVID-19) Envelope Antibody



Species/Host

Rabbit

Reactivity

Virus

Conjugation

Biotin

Tested


ELISA, IHC

Applications

Immunohistochemistry
Validation of SARS-...

Immunogen

Anti-SARS-CoV-2 (COVID-19, 2019-nCoV) Envelope antibody (3531) was raised against a peptide corresponding to 10 amino acids near the amino terminus of SARS-CoV-2 (COVID-19, 2019-nCoV) Envelope protein. The immunogen is located within the first 50 amino acids of SARS-CoV-2 (COVID-19, 2019-nCoV) Envelope.



ELISA Validation.
Coating Antigen:
immun...

Target

E

Preservatives

SARS-CoV-2 (COVID-19, 2019-nCoV) Envelope Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

SARS-CoV-2 COVID-19, 2019-nCoV) Envelope antibody can be stored at 4 ° C for three months and -20 ° C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

IHC: 2 µg/mL Antibody validated: Immunohistochemistry in COVID-19 patient samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID

QHD43418

NCBI

QHD43418

Dilution Range

IHC: 2 µg/mL Antibody validated: Immunohistochemistry in COVID-19 patient