

Product Datasheet

SARS-CoV-2 (COVID-19) Spike Antibody (HRP) (orb1239044)

Description

SARS-CoV-2 (COVID-19) Spike Antibody (HRP)

COVID-19 (+) COVID-19 (-)

Species/Host

Rabbit

Reactivity

Virus

Conjugation

HRP

Tested

ELISA, IHC

Applications

Immunohistochemistry
Validation of SARS-...

Immunogen

Anti-SARS-CoV-2 (COVID-19) Spike antibody (HRP) (3525) was raised against a peptide corresponding to 20 amino acids near the carboxy terminus of SARS-CoV-2 (COVID-19) Spike glycoprotein. The immunogen is located within the last 50 amino acids of SARS-CoV-2 (COVID-19) Spike protein.

Target

S

Preservatives

SARS-CoV-2 (COVID-19) Spike Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

SARS-CoV-2 (COVID-19) Spike antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

IHC: 1 µg/mL Antibody validated:
Immunohistochemistry in COVID-19 patient samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID



[P0DTC2](#)

NCBI

[QHD43416](#)

Dilution Range

IHC: 1 µg/mL Antibody validated:
Immunohistochemistry in COVID-19 patient samples. All other applications and species not yet tested.



ELISA Validation.
Coating Antigen:
immun...