

Product Datasheet

**SARS-CoV-2 (COVID-19) Spike 156-157EFdel
Antibody (Delta Variant) (orb1239050)**

Description

SARS-CoV-2 (COVID-19) Spike 156-

Species/Host

Rabbit

Reactivity

Virus

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Anti-SARS-CoV-2 Spike 156-157EFdel (Delta Variant) antibody (orb1239050) was raised against a peptide corresponding to 11 amino acids near the amino terminal domain of SARS-CoV-2 Delta variant (B.1.617.2) Spike protein. The immunogen is located within 130-180 amino acids of SARS-CoV-2 Delta variant (B.1.617.2) Spike protein.

Target

S

Preservatives

SARS-CoV-2 Spike Delta Variant antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

SARS-CoV-2 Spike Delta Variant antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID


[P0DTC2](#)

NCBI


[QHD43416](#)

Dilution Range


WB: 1 µg/mL; IHC: 0.5 µg/mL.
Antibody validated: Western Blot and IHC in human samples. Anti-SARS-CoV-2 Spike 156-157EFdel (Delta



SARS-Cov-2 Spike
156-157EFdel (Delta
Var...



Immunohistochemistry
Validation of Spike...



WB Validation of of
Delta Variant Spike ...