




Product Datasheet

ORF9c Antibody (orb1239056)

Description	ORF9c Antibody	
Species/Host	Rabbit	 Western Blot Validation with SARS-CoV-2 ...
Reactivity	Virus	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Anti-SARS-CoV-2 (COVID-19) ORF9c antibody (orb1239056) was raised against a peptide corresponding to 16 amino acids near the amino-terminus of SARS-CoV-2 (COVID-19) ORF9c protein. The immunogen is located within the first 50 amino acids of the SARS-CoV-2 (COVID-19) ORF9c protein.	
Target	ORF9c	 Western Blot Validation with SARS-CoV-2 ...
Preservatives	SARS-CoV-2 (COVID-19) ORF9c antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	SARS-CoV-2 (COVID-19) ORF9c antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	 ELISA Validation. Antibodies: SARS-CoV-2...
Isotype	IgG	
Clonality	Polyclonal	
MW	8kD	
Uniprot ID	P0DTD3	
NCBI	P0DTD3	
Dilution Range	WB: 0.25-1 µg/mL Antibody validated: SARS-CoV-2 (COVID-19) ORF9c antibody can detect 2 ng of free peptide at 1 µg/mL in ELISA. It can detect SARS-CoV-2 ORF9c in transfected human cells by WB. All other applications and species not yet tested	
Expiration Date	12 months from date of receipt.	