



Product Datasheet

ORF8 Antibody (orb1239057)

Description	ORF8 Antibody	KDa 1 2
Species/Host	Rabbit	
Reactivity	Virus	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Anti-SARS-CoV-2 (COVID-19) ORF8 antibody (orb1239057) was raised against a peptide corresponding to 13 amino acids near the central domain of SARS-CoV-2 (COVID-19) ORF8 protein. The immunogen is located within 30-80 amino acids of the SARS-CoV-2 (COVID-19) ORF8 protein.	Western Blot Validation with SARS-CoV-2 ...
Target	ORF8	
Preservatives	SARS-CoV-2 (COVID-19) ORF8 antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	SARS-CoV-2 (COVID-19) ORF8 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	ELISA Validation. Antibodies: SARS-CoV-2...
Note	For research use only	
Isotype	IgG	
Clonality	Polyclonal	
MW	14kD	
Uniprot ID	P0DTC8	
NCBI	P0DTC8	
Dilution Range	WB: 1 µg/mL Antibody validated: SARS-CoV-2 (COVID-19) ORF8 antibody can detect 2 ng of free peptide at 1 µg/mL in ELISA. It can detect SARS-CoV-2 ORF8 recombinant protein by WB. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	