

Product Datasheet

Non-structural protein 9 Antibody (orb1239068)

Description

Non-structural protein 9 Antibody

Species/Host

Rabbit

Reactivity

Virus

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Anti-SARS-CoV-2 (COVID-19) NSP9 antibody (orb1239068) was raised against a peptide corresponding to 14 amino acids near the carboxyl terminus of SARS-CoV-2 (COVID-19) NSP9 protein. The immunogen is located within the last 50 amino acids of the SARS-CoV-2 (COVID-19) NSP9 protein

Target

Non-structural protein 9

Preservatives

SARS-CoV-2 (COVID-19) NSP9 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

SARS-CoV-2 (COVID-19) NSP9 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID

[P0DTC1](#)

NCBI


[YP_009742616.1](#)

Dilution Range


WB: 0.5 µg/mL Antibody validated: SARS-CoV-2 (COVID-19) NSP9 antibody can detect 2 ng of free peptide at 1 µg/mL in ELISA.

Expiration Date

12 months from date of receipt.



ELISA
Validation
with SARS-
CoV-2
(COVID-...



Western
Blot
Validation
with SARS-
CoV-2 ...