

Product Datasheet

HEL Antibody (orb1239078)

Description

HEL Antibody

Species/Host

Rabbit

Reactivity

Virus

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications
Immunogen

Anti-SARS-CoV-2 (COVID-19) NSP13 (Helicase)antibody (orb1239078) was raised against a peptide corresponding to 14 amino acids near the center of the SARS-CoV-2 (COVID-19) NSP13 (Helicase) protein. The immunogen is located between 440-490 amino acids of SARS-CoV-2 (COVID-19) NSP13 (Helicase).

Target

HEL

Preservatives

SARS-CoV-2 (COVID-19) NSP13 (Helicase) Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

SARS-CoV-2 COVID-19) NSP13 (Helicase) antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality




Polyclonal

Uniprot ID
P0DTC1
NCBI
YP_009725308.1
Dilution Range

IHC: 0.5 µg/mL. Antibody validated: Immunohistochemistry in human samples. SARS-CoV-2 (COVID-19) NSP13 antibody can detect 2 ng of free peptide at 1 µg/mL in ELISA. It can detect SARS-CoV-2 NSP13 recombinant protein by ELISA. All other applications and species not yet tested

Expiration Date

12 months from date of receipt.


 Immunohistochemistry
Validation of SARS-...

 ELISA Validation with
SARS-CoV-2 (COVID-...

 ELISA Validation.
Antibodies: SARS-CoV-
2...