

Product Datasheet

SARS-CoV-2 (COVID-19) Membrane Antibody (HRP) (orb1239083)

Description

SARS-CoV-2 (COVID-19) Membrane Antibody

COVID-19 (+) COVID-19 (-)

Species/Host

Rabbit

Reactivity

Virus

Conjugation

HRP

Tested


ELISA, IHC

Applications

Immunohistochemistry
Validation of SARS...

Immunogen

Anti-SARS-CoV-2 (COVID-19) Membrane antibody (HRP) (9157) was raised against a peptide corresponding to 13 amino acids near the center of SARS-CoV-2 (COVID-19) Membrane protein. The immunogen is located between 130-180 amino acids of the SARS-CoV-2 (COVID-19) Membrane protein.



Target

M

ELISA

Preservatives

SARS-CoV-2 (COVID-19) Membrane Antibody is supplied in PBS containing 0.02% sodium azide.

Validation.Coating
Antigen: immuno...

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

SARS-CoV-2 (COVID-19) Membrane antibody can be stored at 4 ° C for three months and -20 ° C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

IHC: 2 µg/mL Antibody validated:
Immunohistochemistry in COVID-19 patient samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID

P0DTC5

NCBI

QJW69073

Dilution Range

IHC: 2 µg/mL Antibody validated:
Immunohistochemistry in COVID-19 patient samples. All other applications and species