

Product Datasheet

ARID1A Antibody (orb1239088)

Description
ARID1A Antibody
Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications
Immunogen

ARID1A antibody was raised against an 18 amino acid peptide near the amino terminus of human ARID1A. The immunogen is located within amino acids 580 - 630 of ARID1A.

Target

ARID1A

Preservatives

ARID1A antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

ARID1A antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

ARID1A antibody can be used for detection of ARID1A by Western blot at 1 - 2 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in mouse samples and Immunofluorescence in mouse samples. All other applications and species not yet tested.

Isotype

IgG

Clonality


Polyclonal

MW


Predicted: 133, 227, 251 kDa Observed: 130, 260 kDa

Uniprot ID
[O14497](https://www.uniprot.org/entry/O14497)
NCBI
[NP_006006](https://www.ncbi.nlm.nih.gov/nuccore/NP_006006)
Dilution Range


ARID1A antibody can be used for detection of ARID1A by Western blot at 1 - 2 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in mouse samples and Immunofluorescence in mouse samples. All other applications and species not yet tested.

A B


Western blot analysis of ARID1A in SK-N-BE...



Immunohistochemistry of ARID1A in mouse brain...



Immunofluorescence of ARID1A in mouse brain...