



Product Datasheet

ARPC1A antibody (orb12391)


Description	Goat polyclonal to ARPC1A.
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ARPC1A
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoforms NP_006400.2;NP_001177925.1. No cross reactivity expected with ARPC1B.
Protein Sequence	DFKTLESSIQGLR
Clonality	Polyclonal
MW	41.6; 39.7
Entrez	10552
Dilution Range	ELISA: 1:64000, WB: 0.1-0.3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa lysate using ARPC1A antibody.



Western blot analysis of fetal mouse brain lysate using ARPC1A antibody.



Immunohistochemical staining of human HeLa cells using ARPC1A antibody.