




Product Datasheet

AIFM1 Antibody (orb1239168)

Description	AIFM1 Antibody	
Species/Host	Rabbit	Western Blot Validation in Human Cell Li...
Conjugation	Unconjugated	
Tested Applications	ELISA, ICC, WB	
Immunogen	Anti-AIF antibody (orb1239168) was raised against a peptide corresponding to 14 amino acids near the carboxy terminus of human AIF. The immunogen is located within the last 50 amino acids of AIF.	
Target	AIFM1	
Preservatives	AIF Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/ml	
Storage	AIF antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	Independent Antibody Validation (IAV) vi...
Note	For research use only	
Isotype	IgG	
Clonality	Polyclonal	
MW	Predicted: 67kD Observed: 68kD	
Uniprot ID	O95381	Western Blot Validation in Mouse and Rat...
NCBI	O95381	
Dilution Range	WB: 0.5-2 µg/mL; ICC: 5 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples; Immunocytochemistry in human samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	