

Product Datasheet

PJA1 Antibody (orb1239218)

Description

PJA1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, IF, IHC-P, WB

Immunogen

PJA1 antibody was raised against a 19 amino acid peptide near the amino terminus of human PJA1. The immunogen is located within amino acids 50 - 100 of PJA1.

Target

PJA1

Preservatives

PJA1 antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

PJA1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

PJA1 antibody can be used for detection of PJA1 by Western blot at 0.25 - 0.5 µg/ml. Antibody can also be used for immunohistochemistry starting at 5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in mouse samples and Immunofluorescence in mouse samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 71 kDa Observed: 71 kDa

Uniprot ID


[Q8NG27](#)

NCBI


[NP_660095](#)

Dilution Range


PJA1 antibody can be used for detection of PJA1 by Western blot at 0.25 - 0.5 µg/ml. Antibody can also be used for immunohistochemistry starting at 5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in mouse samples and Immunofluorescence in mouse samples. All other applications and species not yet tested.



Western blot analysis of PJA1 in human brain tissue.



Immunohistochemistry of PJA1 in mouse brain tissue.



Immunofluorescence of PJA1 in mouse brain tissue.