

Product Datasheet

PHF14 Antibody (orb1239241)

Description

PHF14 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, WB

Applications

Immunogen

PHF14 antibody was raised against a 19 amino acid peptide near the amino terminus of human PHF14. The immunogen is located within amino acids 210 - 260 of PHF14 .

Target

PHF14

Preservatives

PHF14 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

PHF14 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

PHF14 antibody can be used for detection of PHF14 by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 98 kDa Observed: 100kDa

Uniprot ID

[O94880](#)

NCBI


[NP_055475](#)

Dilution Range


PHF14 antibody can be used for detection of PHF14 by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Western blot analysis of PHF14 in 293 cells...



Immunofluorescence of PHF14 in 293 cells...