

Product Datasheet

TNFSF13B Antibody (orb1239254)

Description

TNFSF13B Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, ICC, IF, IHC, WB

Immunogen

Anti-BAFF antibody (orb1239254) was raised against a peptide corresponding to 16 amino acids near the carboxy terminus of human BAFF. The immunogen is located within the last 50 amino acids of BAFF.

Target

TNFSF13B

Preservatives

BAFF Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

BAFF antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 31kD Observed: 35 kD

Uniprot ID

[Q9Y275](#)

NCBI


[NP_006564](#)

Dilution Range


WB: 0.25-1 µg/mL; ICC: 1 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples; Immunocytochemistry in human samples; Immunofluorescence in human samples; Immunohistochemistry in human and mouse samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Western Blot Validation in Human HL60 Ce...



Western Blot Validation in Human, Mouse ...



Western Blot Validation with Recombinant...