

Product Datasheet

BCL2L11 Antibody (orb1239259)

Description

BCL2L11 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

Anti-BIM antibody (orb1239259) was raised against a peptide corresponding to 20 amino acids near the amino terminus of human BIM. The immunogen is located within the first 50 amino acids of BIM.

Target

BCL2L11

Preservatives

Bim Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/ml

Storage

Bim antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 22kD Observed: 22kD

Uniprot ID

043521

NCBI


043521

Dilution Range


WB: 0.5 µg/mL; IHC-P: 20 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples; Immunohistochemistry in human samples; Immunofluorescence in human samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Western Blot Validation in Human Cell Li...



Independent Antibody Validation (IAV) vi...



Western Blot Validation in Human Tissue....