

Product Datasheet

CASP9 Antibody (orb1239265)

Description

CASP9 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, IP, WB

Applications

Immunogen

Anti-Caspase 9 antibody (orb1239265) was raised against a peptide corresponding to 20 amino acids near the center of human Caspase 9. The immunogen is located within amino acids 290-340 of Caspase 9.

Target

CASP9

Preservatives

Caspase-9 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Caspase-9 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 46kD
Observed: 46 kD

Uniprot ID

[P55211](#)

NCBI


[P55211](#)

Dilution Range


WB: 1 µg/mL; ICC: 2- 5 µg/mL; IF: 5-20 µg/mL. Antibody validated: Western Blot in human and mouse samples, Immunocytochemistry, Immunofluorescence in human samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Independent
Antibody
Validation
(IAV) vi...



Western Blot
Validation in
Human Cell
Li...



Western Blot
Validation in
Mouse Cell
Li...