


Product Datasheet

PDGFB Antibody (orb1239269)


Description

PDGFB Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Rabbit polyclonal PDGF-B antibody was raised against a 19 amino acid peptide near the carboxy terminus of human PDGF-B. The immunogen is located within the last 50 amino acids of PDGF-B.
Target	PDGFB
Preservatives	PDGF-B Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	PDGF-B antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Note	For research use only
Application notes	PDGF-B antibody can be used for detection of PDGF-B by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in rat samples and Immunofluorescence in human samples. All other applications and species not yet tested.
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 27 kDa
Uniprot ID	P01127
NCBI	NP_002599
Dilution Range	PDGF-B antibody can be used for detection of PDGF-B by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in rat samples and Immunofluorescence in human samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.



Western blot analysis of PDGF-B in rat l...



Immunofluorescence of PDGF-B in human li...