




Product Datasheet

BIRC2 Antibody (orb1239295)

Description	BIRC2 Antibody	
Species/Host	Rabbit	 Western Blot Validation in Human Lung Ly...
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, IF, IHC-P, WB	
Immunogen	Anti-cIAP antibody (orb1239295) was raised against a peptide corresponding to 14 amino acids near the carboxy terminus of human cIAP. c-IAP antibody detects both c-IAP1 and c-IAP2. The immunogen is located within amino acids 540-590 of cIAP.	
Target	BIRC2	 Independent Antibody Validation (IAV) vi...
Preservatives	cIAP Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	cIAP antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	 Western Blot Validation in Human Cell Li...
Isotype	IgG	
Clonality	Polyclonal	
MW	Predicted: 70kD Observed: 70kD	
Uniprot ID	Q13490	
NCBI	NP_001157	
Dilution Range	WB: 1-4 µg/mL; IHC: 10 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples; Immunohistochemistry in human samples; Immunofluorescence in human samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	