


Product Datasheet

C1QTNF3 Antibody (orb1239313)


Description

C1QTNF3 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	Anti-CTRP3 antibody (orb1239313) was raised against a peptide corresponding to 16 amino acids near the center of human CTRP3. The immunogen is located within amino acids 90-140 of CTRP3.
Target	C1QTNF3
Preservatives	CTRP3 Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	CTRP3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 27kD Observed: 30 kD (Post-modification: 1 N-linked glycosylation)
Uniprot ID	Q9BXJ4
NCBI	NP_852100
Dilution Range	WB: 1 µg/mL; IHC: 2 µg/mL; IF: 10 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples; Immunohistochemistry in mouse samples; Immunofluorescence in mouse samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.



Western Blot Validation in Human and Mou...



Western Blot Validation with Recombinant...



Western Blot Validation in Mouse and Rat...