

Product Datasheet

OLIG1 Antibody (orb1239320)

Description

OLIG1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

ELISA, IF, IHC-P, WB

Immunogen

OLIG1 antibody was raised against an 18 amino acid peptide near the amino terminus of human OLIG1. The immunogen is located within amino acids 50 - 100 of OLIG1.

Target

OLIG1

Preservatives

OLIG1 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

OLIG1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

OLIG1 antibody can be used for detection of OLIG1 by Western blot at 1 - 2 µg/ml. Antibody can also be used for Immunohistochemistry at 5 µg/mL. For Immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 30 kDa Observed: 29 kDa

Uniprot ID

[Q8TAK6](#)

NCBI


[NP_620450](#)

Dilution Range


OLIG1 antibody can be used for detection of OLIG1 by Western blot at 1 - 2 µg/ml. Antibody can also be used for Immunohistochemistry at 5 µg/mL. For Immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

31-
24-
14-
7-
← OLIG1

Western blot analysis of OLIG1 in A549 c...



Immunohistochemistry of OLIG1 in human b...



Immunofluorescence of OLIG1 in human bra...