

Product Datasheet

TNFRSF25 Antibody (orb1239336)

Description

TNFRSF25 Antibody

Species/Host

Rabbit

Conjugation

Unconjugated

Tested Applications

ELISA, ICC, IF, WB

Immunogen

DR3 antibody was raised against a peptide corresponding to amino acids in extracellular domain of human DR3 precursor. The immunogen is located within amino acids 50 - 100 of DR3.

Target

TNFRSF25

Preservatives

DR3 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/ml

Storage

DR3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

DR3 antibody can be used for detection of DR3 expression by Western blot at 1 µg/mL. 59 kDa band should be detected. Antibody can also be used for immunocytochemistry starting at 10 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunocytochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

59 kDa

Uniprot ID


[Q93038](#)

NCBI


[AAQ88676](#)

Dilution Range


DR3 antibody can be used for detection of DR3 expression by Western blot at 1 µg/mL. 59 kDa band should be detected. Antibody can also be used for immunocytochemistry starting at 10 µg/mL. For immunofluorescence start at 20 µg/mL.



Western blot analysis of DR3 in Jurkat t...



Immunocytochemistry of DR3 in Jurkat cel...



Immunofluorescence of DR3 in Jurkat cell...