

Product Datasheet

TNFRSF10A Antibody (orb1239340)

Description

TNFRSF10A Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, IHC-P, WB

Applications

Immunogen

Anti-DR4 antibody (orb1239340) was raised against a peptide corresponding to 19 amino acids near the carboxy terminus of human DR4. The immunogen is located within the last 50 amino acids of DR4.

Target

TNFRSF10A

Preservatives

DR4 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

DR4 antibody can be stored at 4°C up to one year. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 50kD Observed: 55kD (Post-modification: 1 N-linked glycosylation)

Uniprot ID

O00220

NCBI


AAC51226

Dilution Range


WB: 1 µg/mL; IHC-P: 10 µg/mL; IHC: 10 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples; Immunohistochemistry in human samples; Immunocytochemistry in human samples and Immunofluorescence in human and rat samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



KO Validation in HeLa Cells. Loading: 10...



Independent Antibody Validation (IAV) vi...



Western Blot Validation in Cell Lines. L...