

Product Datasheet


MKRN3 Antibody (orb1239434)

Description


MKRN3 Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	MKRN3 antibody was raised against a 17 amino acid peptide near the carboxy terminus of human MKRN3. The immunogen is located within amino acids 400 - 450 of MKRN3.
Target	MKRN3
Preservatives	MKRN3 antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	MKRN3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Note	For research use only
Application notes	MKRN3 antibody can be used for detection of MKRN3 by Western blot at 1 - 2 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 55 kDa Observed: 55 kDa
Uniprot ID	Q6NSB6
NCBI	NP_005655
Dilution Range	MKRN3 antibody can be used for detection of MKRN3 by Western blot at 1 - 2 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.

A B



Western blot analysis of MKRN3 in human ...



Immunofluorescence of MKRN3 in human spl...