




Product Datasheet

IKBKE Antibody (orb1239458)

Description	IKBKE Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, IF, IHC-P, WB	
Immunogen	Anti-IKK epsilon antibody (orb1239458) was raised against a peptide corresponding to 16 amino acids near the carboxy terminus of human IKK epsilon. The immunogen is located within the last 50 amino acids of IKK epsilon.	
Target	IKBKE	Induction Validation in Mouse Cell Line.... 
Preservatives	IKK epsilon Antibody is supplied in PBS containing 0.02% sodium azide.	Western Blot Validation in Human and Mou... 
Form/Appearance	Liquid	Western Blot Validation in Mouse Tissue ...
Concentration	1 mg/ml	
Storage	IKK epsilon antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Isotype	IgG	
Clonality	Polyclonal	
MW	Predicted: 80kD Observed: 80 kD	
Uniprot ID	Q14164	
NCBI	AF241789	
Dilution Range	WB: 1-2 µg/mL; IHC: 10 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human and mouse samples; Immunohistochemistry in human samples; Immunofluorescence in human samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	