




Product Datasheet

IL17A Antibody (orb1239477)

Description	IL17A Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, IF, IHC-P, WB	
Immunogen	Anti-IL-17 antibody (orb1239477) was raised against a peptide corresponding to 16 amino acids near the center of mature human IL-17. The immunogen is located within amino acids 50-100 of IL-17.	 Western Blot Validation with Recombinant...
Target	IL17A	 Western Blot Validation with Human Splene...
Preservatives	IL-17 Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	IL-17 antibody can be stored at 4°C for up to one year. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Isotype	IgG	
Clonality	Polyclonal	
MW	Predicted: 17.5kDaObserved: 18 kDa	
Uniprot ID	Q16552	
NCBI	NP_002181	
Dilution Range	WB: 0.125-2 µg/mL; IF: 5-20 µg/mL; IHC: 1-5 µg/mL .Antibody validated: Western Blot in human samples; Immunofluorescence in mouse samples; Immunohistochemistry in human and mouse samples. All other applications and species not yet tested.	 Immunohistochemistry Validation of IL-17...
Expiration Date	12 months from date of receipt.	