

Product Datasheet

IL32 Antibody (orb1239494)

Description

IL32 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

Anti-IL-32 antibody (orb1239494) was raised against a peptide corresponding to 15 amino acids near the carboxy terminus of human IL-32. The immunogen is located within the last 50 amino acids of IL-32.

Target

IL32

Preservatives

IL-32 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

IL-32 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 26kD
Observed: 26 kD

Uniprot ID

[P24001](#)

NCBI

[AAH09401](#)


Dilution Range

WB: 0.5-4 µg/mL; IHC: 2-10 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human and mouse samples; Immunofluorescence in human samples. All other applications and species not yet tested.


Expiration Date

12 months from date of receipt.


KDa 1 2 3 4



Western Blot
Validation with
Recombinant...



Western Blot
Validation in Human
Spleen ...



Immunohistochemistry
Validation of IL-32...