


Product Datasheet

IL33 Antibody (orb1239503)

Description


IL33 Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	Anti-IL-33 antibody (orb1239503) was raised against a 17 amino acid synthetic peptide from near the center of human IL-33. The immunogen is located within the 150-200 amino acids of IL-33.
Target	IL33
Preservatives	IL-33 Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	IL-33 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 31kD Observed: 33kD
Uniprot ID	O95760
NCBI	NP_254274
Dilution Range	WB: 1-2 µg/mL; IF: 20 µg/mL; IHC-P: 2.5 µg/mL. Antibody validated: Western Blot in human and mouse samples; Immunofluorescence and Immunohistochemistry in human samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.




Western Blot
Validation of
IL-33 in A20

...



Western Blot
Validation of
IL-33 with
Hu...



Western Blot
Validation of
IL-33 in
Huma...