

Product Datasheet

IL9 Antibody (orb1239508)

Description

IL9 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

ELISA, IF, IHC-P, WB

Immunogen

Anti-IL-9 antibody (orb1239508) was raised against a peptide corresponding to 16 amino acids near the carboxy terminus of human IL-9. The immunogen is located within the last 50 amino acids of IL-9.

Target

IL9

Preservatives

IL-9 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

IL-9 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 16kD
Observed: 22 kD (Post-modification: 4 N-linked glycosylations)

Uniprot ID

[P15248](#)

NCBI


[NP_000581](#)

Dilution Range


WB: 1-4 µg/mL; IHC: 5 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples; Immunofluorescence in human samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Western Blot Validation in Human Spleen ...



Western Blot Validation with Recombinant...



Western Blot Validation in Tissue and Ce...