

Product Datasheet

IRF7 Antibody (orb1239524)

Description

IRF7 Antibody

Species/Host

Rabbit

Conjugation

Unconjugated

Tested Applications

ELISA, IHC-P, WB

Immunogen

Anti-IRF7 antibody (orb1239524) was raised against a peptide corresponding to 17 amino acids near the center of human IRF7. The immunogen is located within amino acids 200-230 of IRF7.

Target

IRF7

Preservatives

IRF7 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

IRF7 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID

[Q92985](#)

NCBI


[Q92985](#)

Dilution Range


WB: 0.5-1 µg/mL; IF: 20 µg/mL; IHC: 2-5 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples. Immunofluorescence in mouse samples. Immunohistochemistry in human, mouse and rat samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



KO Validation
in A549 Cell
Lysate.
Loadi...



Western Blot
Validation in
Human Cell
Li...



Western Blot
Validation in
Human
Tissues...