

Product Datasheet

KDM3A Antibody (orb1239525)

Description

KDM3A Antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

Anti-JMJD1A antibody (orb1239525) was raised against a peptide corresponding to 16 amino acids near the amino terminus of human JMJD1A. The immunogen is located within amino acids 50-100 of JMJD1A.

Target

KDM3A

Preservatives

JMJD1A Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

JMJD1A antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 147kD Observed: 147kD

Uniprot ID

[Q9Y4C1](#)

NCBI


[NP_060903](#)

Dilution Range


WB: 1-2 µg/mL; IHC: 5 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in rat samples; Immunofluorescence in rat samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Western Blot Validation in Human Cell Li...



Independent Antibody Validation (IAV) vi...



Western Blot Validation in Human HepG2 C...