

Product Datasheet

KL Antibody (orb1239544)

Western blot analysis of the 116 kDa band in various tissues. The blot shows a single band at 116 kDa in all lanes, indicating the presence of the protein in all tissues. The lanes are labeled: 293, A43, A54, CaC, HeF, K56, MCF, Jur, SK-, THF, NIH, and F3. Molecular weight markers are indicated on the left at 116, 116, and 42 kDa.

SDS-PAGE gel showing a single band at approximately 100 kDa, indicated by an arrow. Molecular weight markers are labeled on the left: 250, 130, 95, 72, and 55 kDa.

Western Blot Validation in Rat Kidney Tissue