




Product Datasheet

LIF Antibody (orb1239553)

Description	LIF Antibody	
Species/Host	Rabbit	 Western Blot Validation in Mouse 3T3 Cel...
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, IF, WB	
Immunogen	Anti-LIF antibody (orb1239553) was raised against a peptide corresponding to 16 amino acids near the center of human LIF. The immunogen is located within amino acids 50-100 of LIF.	
Target	LIF	 Western Blot Validation in Human Cell Li...
Preservatives	LIF Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	LIF antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	 Western Blot Validation with Recombinant...
Isotype	IgG	
Clonality	Polyclonal	
MW	Predicted: 22kD Observed: 26kD (Post-modifications: N-linked Glycosylations)	
Uniprot ID	P15018	
NCBI	NP_002300	
Dilution Range	WB: 1-4 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples; Immunofluorescence in mouse samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	