

Product Datasheet

KREMEN1 Antibody (orb1239559)

Description

KREMEN1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

Rabbit polyclonal KREMEN1 antibody was raised against an 18 amino acid peptide near the carboxy terminus of human KREMEN1. The immunogen is located within the last 50 amino acids of KREMEN1.

Target

KREMEN1

Preservatives

KREMEN1 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

KREMEN1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

KREMEN1 antibody can be used for detection of KREMEN1 by Western blot at 0.125 - 0.25 µg/mL. Antibody validated: Western Blot in rat samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

54 kDa

Uniprot ID

Q96MU8

NCBI


NP_114434

Dilution Range


KREMEN1 antibody can be used for detection of KREMEN1 by Western blot at 0.125 - 0.25 µg/mL. Antibody validated: Western Blot in rat samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Western blot analysis of KREMEN1 in rat ...



Immunohistochemistry of KREMEN1 in human...



Immunofluorescence of KREMEN1 in human S...