

Product Datasheet

KHDC1 Antibody (orb1239581)

Description

KHDC1 Antibody

A B

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications

Western blot analysis
of KHDC1 in rat li...

Immunogen

KHDC1 antibody was raised against a 15 amino acid peptide near the amino terminus of human KHDC1. The immunogen is located within amino acids 30 - 80 of KHDC1.

Target

KHDC1

Preservatives

KHDC1 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

KHDC1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

KHDC1 antibody can be used for detection of KHDC1 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in rat samples and Immunohistochemistry in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 18 kDa Observed: 26 kDa

Uniprot ID

[Q4VXA5](#)

NCBI



[NP_085045](#)

Dilution Range

KHDC1 antibody can be used for detection of KHDC1 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in rat samples and Immunohistochemistry in human samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Immunohistochemistry
of KHDC1 in human l...