

Product Datasheet

KANK3 Antibody (orb1239590)

Description

KANK3 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications

Immunogen

KANK3 antibody was raised against an 18 amino acid peptide near the amino terminus of human KANK3 .The immunogen is located within amino acids 50 - 100 of KANK3.

Target

KANK3

Preservatives

KANK3 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

KANK3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

KANK3 Antibody can be used for detection of KANK3 by Western blot starting at 1 µg/mL.Antibody validated: Western Blot in human samples and Immunohistochemistry in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 90 kDa Observed: 100 kDa

Uniprot ID

[Q6NY19](#)

NCBI


[NP_940873](#)

Dilution Range


KANK3 Antibody can be used for detection of KANK3 by Western blot starting at 1 µg/mL.Antibody validated: Western Blot in human samples and Immunohistochemistry in human samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Western blot analysis of KANK3 in A431 c...



Immunohistochemistry of KANK3 in human k...