



## **Product Datasheet**

JPH3 Antibody (orb1239596)



## Descriptionnts. JPH3 Antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

**Conjugation** Unconjugated

Tested ELISA, IF, IHC-P, WB

**Applications** 

Immunogen JPH3 antibody was raised against a 18 amino acid

synthetic peptide near the carboxy terminus of human JPH3. The immunogen is located within amino

acids 540 - 590 of JPH3.

Target JPH3

**Preservatives** JPH3 Antibody is supplied in PBS containing 0.02%

sodium azide.

Form/Appearance Liquid

Concentration 1 mg/mL

**Storage** JPH3 antibody can be stored at 4°C for three months

and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed

to prolonged high temperatures.

**Note** For research use only

**Application notes** JPH3 antibody can be used for detection of JPH3 by

Western blot at 1  $\mu$ g/mL. Antibody can also be used for immunohistochemistry starting at 2.5  $\mu$ g/mL. For immunofluorescence start at 20  $\mu$ g/mL.Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other

applications and species not yet tested.

**Isotype** IgG

**Clonality** Polyclonal


Uniprot ID Q8WXH2

NCBI NP 065706


**Dilution Range** JPH3 antibody can be used for detection of JPH3 by

Western blot at 1  $\mu$ g/mL. Antibody can also be used for immunohistochemistry starting at 2.5  $\mu$ g/mL. For immunofluorescence start at 20  $\mu$ g/mL.Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and


www.biorbyt.com



Western blot analysis of JPH3 in Daudi c...



Immunohistochemistry of JPH3 in human br...



Immunofluorescence of JPH3 in Human Brai...

68 TW Alexander Drive<br/>br>Research Triangle Park<br/>br>Durham, North