

Product Datasheet

JAG1 Antibody (orb1239611)

Description
JAG1 Antibody
Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications
Immunogen

JAG1 antibody was raised against a 19 amino acid peptide near the carboxy terminus of human JAG1. The immunogen is located within the last 50 amino acids of JAG1.

Target

JAG1

Preservatives

JAG1 antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

JAG1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

JAG1 antibody can be used for detection of JAG1 by Western blot at 1 - 2 µg/ml. Antibody can also be used for immunohistochemistry starting at 5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality


Polyclonal

MW

Predicted: 134 kDa Observed: 130 kDa

Uniprot ID
P78504
NCBI
NP_000205
Dilution Range


JAG1 antibody can be used for detection of JAG1 by Western blot at 1 - 2 µg/ml. Antibody can also be used for immunohistochemistry starting at 5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.



Western blot analysis of JAG1 in HeLa ce...



Immunohistochemistry of JAG1 in human sp...



Immunofluorescence of JAG1 in human sple...