




Product Datasheet

MYD88 Antibody (orb1239635)

Description	MYD88 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, IF, IP, WB	
Immunogen	Anti-MYD88 antibody (orb1239635) was raised against a peptide corresponding to 17 amino acids near carboxy terminus of human MYD88 isoform 1. The immunogen is located within the last 50 amino acids of MYD88.	
Target	MYD88	KO Validation in HeLa Cells. Loading: 10... 
Preservatives	MYD88 Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	MYD88 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	Western Blot Validation of MyD88 in huma... 
Isotype	IgG	
Clonality	Polyclonal	
MW	Predicted: 35kD Observed: 35kD	
Uniprot ID	U70451	
NCBI	AAB49967.1	Independent Antibody Validation (IAV) vi...
Dilution Range	WB: 2 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples; Immunofluorescence and Immunoprecipitation in human samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	