

Product Datasheet

MYD88 Antibody (orb1239647)

MYD88 Antibody

RabbitHuman, Mouse, RatUnconjugated

ELISA, IF, IHC-P, IP, WB

Anti-MYD88 antibody (orb1239647) was raised against a peptide corresponding to 16 amino acids near the center of human MYD88 isoform 1. The immunogen is located within amino acids 220 - 270 of MYD88.

MYD88

MYD88 Antibody is supplied in PBS containing 0.02% sodium azide.

Liquid

1 mg/mL

MYD88 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

For research use only

IgG

Polyclonal

Predicted: 35kD Observed: 35kD

U70451AAB49967.1

WB: 1 - 2 µg/mL; IHC-P: 2 µg/mL; IF: 10-20 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples; Immunohistochemistry, Immunoprecipitation and Immunofluorescence in human samples. All other applications and species not yet tested.

12 months from date of receipt.

