

Product Datasheet

IL36A Antibody (orb1239652)

Description

IL36A Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

IL-36A antibody was raised against an 11 amino acid peptide near the carboxy terminus of human IL-36A. The immunogen is located within the last 50 amino acids of IL-36A.

Target

IL36A

Preservatives

IL-36A antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

IL-36A antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

IL-36A antibody can be used for detection of IL-36A by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 17 kDa Observed: 16 kDa

Uniprot ID

[Q9UHA7](#)

NCBI



[NP_055255](#)

Dilution Range


IL-36A antibody can be used for detection of IL-36A by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Immunofluorescence of IL-36A in human br...



Immunohistochemistry of IL-36A in human ...



Western blot analysis of IL-36A in MCF7 ...