




Product Datasheet

NPHS1 Antibody (orb1239658)

Description	NPHS1 Antibody	KDa 1 2
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, IF, IHC-P, WB	
Immunogen	Anti-Nephrin antibody (orb1239658) was raised against a peptide corresponding to 15 amino acids near the center of human Nephrin. The immunogen is located within amino acids 1100-1150 of Nephrin.	
Target	NPHS1	Western Blot Validation in HeLa Cell Lys... 
Preservatives	Nephrin antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	Nephrin antibody can be stored at 4°C for three months and -20°C, stable for up to one year.	
Note	For research use only	Western Blot Validation in Human and Mou... 
Isotype	IgG	
Clonality	Polyclonal	
MW	Predicted: 125, 137 kDa Observed: 125 kDa	
Uniprot ID	O60500	
NCBI	NP_004637	Immunohistochemistry Validation of Nephr...
Dilution Range	WB: 1-2 µg/mL; IHC: 5 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human and mouse samples; Immunohistochemistry in human samples; Immunofluorescence in human samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	