

Product Datasheet

NPTX2 Antibody (orb1239671)

Description

NPTX2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, IF, IHC-P, WB

Immunogen

Anti-NPTX2 antibody (orb1239671) was raised against a peptide corresponding to 16 amino acids near the center of human NPTX2. The immunogen is located within amino acids 170 - 220 of NPTX2.

Target

NPTX2

Preservatives

NPTX2 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

NPTX2 antibody can be stored at 4°C up to one year. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 47kD
Observed: 47 kD

Uniprot ID

[P47972](#)

NCBI


[NP_002514](#)

Dilution Range


WB: 0.5-4 µg/mL; IHC: 5 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples; Immunohistochemistry in human and mouse samples; Immunofluorescence in human and mouse samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Western Blot Validation in Mouse Brain T...



Independent Antibody Validation (IAV) vi...



Western Blot Validation in Mouse Brain T...