

Product Datasheet

STMN1 Antibody (orb1239690)

Description

STMN1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, IF, IHC-P, WB

Immunogen

Anti-op18 antibody (orb1239690) was raised against a peptide corresponding to 17 amino acids near the amino terminus of human op18. The immunogen is located within the first 50 amino acids of op18.

Target

STMN1

Preservatives

op18 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

op18 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 17kD Observed: 17 kD

Uniprot ID

P16949

NCBI


AAH82228

Dilution Range


WB: 0.5-2 µg/mL; IHC: 2.5 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples; Immunohistochemistry and Immunofluorescence in human samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Western Blot Validation in Human Brain T...



Independent Antibody Validation (IAV) vi...



Western Blot Validation in Human Cell Li...