

Product Datasheet

ICOS Antibody (orb1239700)

Description

ICOS Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IF, IHC, WB

Applications

Immunogen

ICOS antibody was raised against a peptide corresponding to 17 amino acids near the carboxy terminus of human ICOS. The immunogen is located within the last 50 amino acids of ICOS.

Target

ICOS

Preservatives

ICOS Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

ICOS antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

ICOS antibody can be used for Immunohistochemistry at 2 µg/mL. For Immunofluorescence start at 20 µg/mL. Antibody validated: Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID

[Q9Y6W8](#)

NCBI



[Q9Y6W8](#)

Dilution Range


ICOS antibody can be used for Immunohistochemistry at 2 µg/mL. For Immunofluorescence start at 20 µg/mL. Antibody validated: Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Immunohistochemistry of ICOS in human lymphoid tissue.



Immunofluorescence of ICOS in human lymphoid tissue.