

Product Datasheet

GRIP1 Antibody (orb1239727)

Description

GRIP1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

Rabbit polyclonal GRIP1 antibody was raised against a 17 amino acid peptide near the carboxy terminus of human GRIP1. The immunogen is located within the last 50 amino acids of GRIP1.

Target

GRIP1

Preservatives

GRIP1 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

GRIP1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

GRIP1 antibody can be used for detection of GRIP1 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in rat samples and Immunofluorescence in rat samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 118 kDa Observed: 115 kDa

Uniprot ID

[Q9Y3R0](#)

NCBI

[NP_066973](#)

Dilution Range

GRIP1 antibody can be used for detection of GRIP1 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in rat samples and Immunofluorescence in rat samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.


A B

130-


95-

◀ GRIP1

Western blot analysis of GRIP1 in 293 ce...



Immunohistochemistry of GRIP1 in rat bra...



Immunofluorescence of GRIP1 in rat brain...