

## Product Datasheet

### GOLGA5 Antibody (orb1239742)

## Description

GOLGA5 Antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse, Rat

### Conjugation

Unconjugated

### Tested

ELISA, IF, IHC-P, WB

### Applications

### Immunogen

GOLGA5 antibody was raised against an 18 amino acid peptide near the carboxy terminus of human GOLGA5. The immunogen is located within amino acids 540 - 590 of GOLGA5.

### Target

GOLGA5

### Preservatives

GOLGA5 Antibody is supplied in PBS containing 0.02% sodium azide.

### Form/Appearance

Liquid

### Concentration

1 mg/mL

### Storage

GOLGA5 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

### Note

For research use only

### Application notes

GOLGA5 antibody can be used for detection of GOLGA5 by Western blot at 1 - 2 µg/mL. Antibody can also be used for Immunohistochemistry starting at 5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

### Isotype

IgG

### Clonality

Polyclonal

### MW

Predicted: 64, 80 kDa Observed: 65, 78 kDa

### Uniprot ID


[Q8TBA6](#)

### NCBI


[NP\\_005104](#)

### Dilution Range


GOLGA5 antibody can be used for detection of GOLGA5 by Western blot at 1 - 2 µg/mL. Antibody can also be used for Immunohistochemistry starting at 5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.



Western blot analysis of GOLGA5 in human...



Immunohistochemistry of GOLGA5 in human...



Immunofluorescence of GOLGA5 in human te...