

Product Datasheet

TNFSF18 Antibody (orb1239752)

Description

TNFSF18 Antibody

Species/Host

Rabbit

Conjugation

Unconjugated

Tested Applications

ELISA, ICC, WB

Immunogen

GITRL antibody was raised against purified recombinant human GITR ligand.

Target

TNFSF18

Preservatives

GITRL Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

GITRL antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

GITRL antibody can be used for the detection of GITRL by Western blot at 1 µg/mL. Antibody can also be used for immunocytochemistry starting at 10 µg/mL. Antibody validated: Immunocytochemistry in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID

[Q9UNG2](#)

NCBI


[Q9UNG2](#)

Dilution Range


GITRL antibody can be used for the detection of GITRL by Western blot at 1 µg/mL. Antibody can also be used for immunocytochemistry starting at 10 µg/mL. Antibody validated: Immunocytochemistry in human samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Western blot analysis of (A) 5 ng, (B) 2...



Immunocytochemistry of GITRL in THP-1 ce...