


Product Datasheet

FGFR3 Antibody (orb1239804)

Description

FGFR3 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	FGFR3 antibody was raised against a 19 amino acid peptide near the carboxy terminus of human FGFR3. The immunogen is located within amino acids 630 - 680 of FGFR3.
Target	FGFR3
Preservatives	FGFR3 Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	FGFR3 antibody can be stored at 4°C up to one year. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	FGF3 antibody can be used for detection of FGFR3 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in rat samples and Immunofluorescence in rat samples. All other applications and species not yet tested.
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 89 kDa Observed: 9+2, 98 kDa
Uniprot ID	P22607
NCBI	NP_001156685
Dilution Range	FGF3 antibody can be used for detection of FGFR3 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in rat samples and Immunofluorescence in rat samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.



Western blot analysis of FGFR3 in SK-N-S...



Immunohistochemistry of FGFR3 in rat bra...



Immunofluorescence of FGFR3 in rat brain...