




Product Datasheet

ANP32A Antibody (orb1239865)

Description	ANP32A Antibody	
Species/Host	Rabbit	
Conjugation	Unconjugated	
Tested Applications	ELISA, ICC, IF, WB	
Immunogen	Anti-PHAP I antibody (orb1239865) was raised against a peptide corresponding to 15 amino acids near the carboxy terminus of human PHAP I. The immunogen is located within the last 50 amino acids of PHAP I.	Western Blot Validation in Human Raji Ce...
Target	ANP32A	
Preservatives	PHAP I Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	PHAP I antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	Independent Antibody Validation (IAV) vi...
Note	For research use only	
Isotype	IgG	
Clonality	Polyclonal	
MW	Predicted: 29kD Observed: 29 kD	
Uniprot ID	P39687	Western Blot Validation in Human Cell Li...
NCBI	P39687	
Dilution Range	WB: 2-4 µg/mL; ICC: 2 µg/mL; IF: 10 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples; Immunocytochemistry in human samples; Immunofluorescence in human samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	