

Product Datasheet

ANP32A Antibody (orb1239868)

Description

ANP32A Antibody

Species/Host

Rabbit

Conjugation

Unconjugated

Tested Applications

ELISA, IF, IHC-P, WB

Immunogen

Anti-PHAP I antibody (orb1239868) was raised against a peptide corresponding to 14 amino acids near the carboxy terminus of human PHAP I . The immunogen is located within the last 50 amino acids of PHAP I.

Target

ANP32A

Preservatives

PHAP I Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

PHAP I antibody can be stored at 4°C up to one year. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 29kD
Observed: 29 kD

Uniprot ID

P39687

NCBI

P39687

Dilution Range

WB: 1 µg/mL; IHC: 2 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples; Immunohistochemistry in mouse samples; Immunofluorescence in mouse samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.

A B C

49

35
28
21

PHAP I

Western Blot Validation in (A) Human Raj...

K562 A431 A549 Daudi HepG2 Hs5787 Jurkat THP-1 U937

Independent Antibody Validation (IAV) vi...

K562 A431 A549 CaCo-2 HepG2 Hs5787 MCF-7 SK-N-SH THP-1 U937 YB2/O

Western Blot Validation in Human, Mouse ...